




BR00150695 | 24 hours | single PFA fixation | 14 cell lines




BR00150696 | 48 hours | single PFA fixation | 14 cell lines



BR00150697 | 72 hours | single PFA fixation | 14 cell lines



BR00150698 | 24 hours | single PFA fixation | 19 cell lines



Seeding density
(cells/well)

- 4,000
- 8,000
- 12,000


Plate coating

- Standard
- Laminin
- Synthemax


No cells

- Media only
- Empty

BR00150699 | 48 hours | single PFA fixation | 19 cell lines



BR00150700 | 72 hours | single PFA fixation | 19 cell lines



Seeding density
(cells/well)

- 1,000
- 4,000
- 8,000
- 12,000


Plate coating

- Standard
- Laminin
- Synthemax

No cells

- Media only
- Empty

BR00150701 | 24 hours | double PFA fixation | 25 cell lines



Seeding density
(cells/well)

- 1,000
- 2,000
- 4,000
- 8,000
- 12,000


Plate coating

- Standard
- Laminin
- Synthemax


No cells

- Media only
- Empty

BR00150702 | 48 hours | double PFA fixation | 25 cell lines



BR00150703 | 72 hours | double PFA fixation | 25 cell lines



Seeding density
(cells/well)

- 1,000
- 2,000
- 4,000
- 8,000
- 12,000

Plate coating

- Standard
- Laminin
- Synthemax

No cells

- Media only
- Empty